



# PARASOL

## Directions for Use

The Parasol Wired Falls Prevention Seatbelt is exclusively used with the Parasol Wired Falls Monitor to alert the staff when a patient is attempting to exit a chair. The red loop allows the patient and staff to easily unlock the belt with a simple sweeping motion that does not require individual finger dexterity. This device is used to alert the nurse staff and **should not** be used as a restraint.

**WARNING:** Before using this device with a patient, read through the following instructions in their entirety. Ensure that the device is used correctly and oriented properly to ensure proper response in the event of the patient attempting to exit their chair. Examine the device for damage before use. **Do not** use the device if there is any damage that effects its functionality. **Always** test the device using the following instructions before using with a patient. Failure to use listed instructions or using a damaged device can cause serious injury to patient. This device is design for use with Parasol Wired Falls Prevention Monitor only.

## Setting Up the Seatbelt

### Attaching Seatbelt to Chair:

Begin by looping the seatbelt around a solid support, recommend using a chair leg or a solid bar, and clip each side onto itself by clipping one of the d-hooks on the green section. See **Figure 1**. Leave the other side unclipped until patient is positioned in the chair. Gently pull on seatbelt to verify that the clip is properly connected and connection location can support the patient's weight.



Figure 1

### Fitting Seatbelt for Patient:

Position the patient in the chair. Guide the unclipped side of the seatbelt between the armrest and the back of the chair to obtain proper securement around patient. Secure the seatbelt around the patient and guide the seatbelt between the arm rest and the back of the chair on the other side. See **Figure 2**. Attach the other metal clip in the the same manner as the first by looping it around a solid support and clipping it onto itself on the green section.



Figure 2

Adjust the seatbelt length by sliding the buckle to shorten or increase the length. See **Figure 3**.




Figure 3

## Connecting Seatbelt to Parasol Wired Falls Monitor

Insert Parasol Wired Falls Monitor into chair mount. Gently pull on monitor to verify properly secured in mount. Insert cord from seatbelt into either the Pad1 or Pad2 outlets on the monitor, See **Figure 4**. Guide the cord so it won't become tangled in the wheel of the wheelchair. The monitor will beep once when plug is inserted into monitor (if seatbelt is engaged).



Figure 4

## Using Seatbelt with Parasol Wired Falls Monitor

To remove the seatbelt, pull the red loop with a simple sweeping motion. See **Figure 5**. The Parasol Wired Falls Prevention Monitor will alarm when seatbelt is disengaged, Velcro blue buckle is separated. See **Figure 6**. To silence the alarm, press the RESET button on the monitor or re-engage the seatbelt. When the seatbelt is re-engaged, the monitor will beep once to confirm connection.



Figure 5



Figure 6

## Cleaning and on-going care of your Seatbelt MicrobeCare™ Treatment on Parasol Monitors and Wireless Pads



For any gross contamination, food, drink etc. or between use with different patients, the seatbelt can be cleaned by wiping it with a wet towelette using water or any standard Hospital disinfectant.

**Continuous Built-in protection** - for your convenience, all Parasol seatbelts are treated with a permanently bonded *antimicrobial* barrier called MicrobeCare™. MicrobeCare™ is non-leaching and has proven viricidal, bactericidal, germicidal and tuberculocidal efficacy. A list of approved cleaning agents and a list of microbes against which the MicrobeCare™ technology has been tested may be obtained by visiting the MicrobeCare™ website at [www.MicrobeCare.com](http://www.MicrobeCare.com).

For a complete falls prevention system, add Parasol Pressure Pads and Incontinence Monitoring System to your Parasol Wired Falls Seatbelt. For a complete list of Parasol Falls Prevention products visit <http://www.parasolmed.com/falls-prevention/> or call Parasol Medical LLC.



Tel: 847-520-3000  
Fax: 847-520-1234  
Email: [info@parasolmed.com](mailto:info@parasolmed.com)

Parasol Medical LLC  
1602 Barclay Blvd.  
Buffalo Grove, IL 60089